

# INFĒRA<sup>®</sup> mini



Portable Light Therapy Device  
USER MANUAL

## Features:

The INFĒRA mini is a compact device that utilizes four different technologies to improve the skin around the delicate eye contour.

**Red Light Therapy** mode that deeply penetrates to stimulate collagen production, improve firmness, and combat discoloration – targeting the appearance of wrinkles and fine lines.

**Blue Light Therapy** mode minimizes swelling and puffiness, combats discoloration, diminishes dark circles and promotes a tighter and more refined surface.

**High-Frequency Vibrations** gently and effectively stimulate dermal layers.

**Thermal Heat** Therapy Assists in stimulating the blood flow and helps boost the absorption of your

## How to Use:

1. Wash face and pat dry.
2. Apply a skincare product of choice, such as a serum, lotion, or cream.
3. Power on device by selection from the following functions:

### Red Light Therapy

- With micro-vibration and heat: Press the left button once (1x).
- With micro-vibration, heat and ion therapy: Press the left button once (2x).
- Without micro-vibration and ion therapy: Press the left button three times (3x).

### Blue Light Therapy

- With micro-vibration: Press the right button once (1x).
- With micro-vibration and ion therapy: Press the right button once (2x).
- Without micro-vibration and ion therapy: Press the right button three times (3x).

4. Once function has been selected, massage the eye area for about 20 to 30 seconds each eye.
5. For best results, use daily. Blue light AM, Red light PM

## How to clean:

1. Ensure that the device is powered off.
  2. Wipe the device with a dry, clean cloth.
- Note: Do not use any detergents or water.
3. Store the device in a dry and ventilated area.

### Cautions:

- Please read the device manual before using and save for future reference.
- Avoid direct open eye contact with light and device applicators.
- Keep out of reach of children.
- For external use only.
- Stop using the device immediately if you feel any signs of pain or discomfort.
- The device is not waterproof. Please do not place it in water or use it on wet skin.
- Do not apply the device onto other sensitive areas of the body.
- This device is not recommended for use by pregnant women or a person wearing a pacemaker.
- Slight redness may occur after immediate use of the device. If redness persists after 24 hrs, discontinue use.

## Conditions for storage:

- Ensure the device is powered off and completely cool before storing.
- Ensure the device is placed in a chemical free environment.
- Ensure the device is out of reach from children.
- Store the device in a dry, clean area.

## Warranty

This device includes a three year warranty period, effective from date of purchase. The warranty scope refers to the issues with INFĒRA products, due to production, technique or quality. In the following circumstances, user is responsible, even during the warranty period:

- Damage due to usage, maintenance or storage not according to user manual
- Damage due to incorrect operation
- Damage due to unauthorized service, not approved by INFĒRA
- Damage due to force majeure (such as lightning, earthquakes, fire, flood, etc.)

You must have a valid purchase receipt and register your warranty on our website after purchase. INFĒRA has the right to void your warranty if you provide any incorrect or altered information.

## Limited Warranty

We offer a limited three year warranty for your INFĒRA mini device. If your product has any manufacturing or workmanship defects, INFĒRA or its official distributors will fix or replace the device.

Warranty will NOT cover normal wear and tear of the device.

Warranty will NOT cover misuse or mistreatment of the device.

Register your new device at:

[www.INFERA-US.com](http://www.INFERA-US.com)

To protect your new INFĒRA against any other accidental damage, please visit our store and purchase the INFĒRA Care.

[www.INFERA-US.com/STORE](http://www.INFERA-US.com/STORE)

If your INFĒRA infrared device is broken, damaged, in need of repair, or any other user assistance, please contact customer service:

[info@infera-us.com](mailto:info@infera-us.com) or call 866-698-2181



# Warranty Card

## ▼ Fill in by customer

Name	_____
Purchase Date	_____
Phone	_____
Address	_____

1. The warranty period is three years, effective from the date of purchase.
2. The warranty scope refers to the issues with INFĒRA products, due to production, technique or quality.
3. In the following circumstances, user is responsible, even during the warranty period:
  - a) Damage due to usage, maintenance or storage not according to user manual.
  - b) Damage due to incorrect operation.
  - c) Damage due to unauthorized service, not approved by INFĒRA
  - d) Warranty card and valid purchase invoice. Must be presented.
  - e) Incorrect or altered record of warranty information.
  - f) Damage due to force majeure (such as lightning, earthquakes, fire, flood, etc.)
4. When contacting brand distributor or manufacturer for warranty repair, please present this warranty card and purchase invoice.
5. During warranty period, the manufacturer will replace the faulty part free.  
Shipping is the responsibility of the customer.
6. The warranty card shall be filled in clear, block letters with no alteration allowed.  
Please keep the warranty card. If lost, there will be no duplicate provided.
7. The final interpretation of the warranty card belongs to INFĒRA.
8. Due to product upgrade or package upgrade, the information contained in this warranty is subject to change without notice.

Contact us: Email: [info@infera-us.com](mailto:info@infera-us.com), Phone: 866-698-2181  
[www.INFERA-US.com](http://www.INFERA-US.com)